

Setting Up CareFlow Check-In on your iPad

Tap **Settings > General > Accessibility > Guided Access** to set up Guided Access. From there you can:

- Turn Guided Access on or off
- Set a passcode that controls the use of Guided Access and prevents someone from leaving an active session
- Set whether the device can go to sleep during a session

Start a Guided Access session

- Open the app you want to run.
- Triple-click the Home button.
- Adjust settings for the session, then click Start.

Disable app controls and areas of the app screen (Image Step 1)

- Circle any part of the screen you want to disable. You should look to disable all browser buttons above the CareFlow Logo on the screen.
- Use the handles to adjust the area.



Keep iPhone, iPad, or iPod touch from switching from portrait to landscape or from responding to any other motions (Image Step 2)

- Turn off Motion.

End a Guided Access session

- Triple-click the Home button.
- Enter the Guided Access passcode.